

DSC-ETS STIRRUP USER INSTRUCTIONS

- 1)**Installation:** Fully extend the leg rests on the exam table you intend to use. Place the wedge included with your set on top or under the footrest. If you place the wedge on top of the footrest it will angle the Stirrup to accommodate shorter legs and if you place the wedge under the footrest it will accommodate longer legs. Either positioning will accommodate most patients but if you have an exceptionally petite or tall patient they can be adjusted to accommodate a wider range.
- 2)Loosen the velcro strap from one side of the Stirrup base.
- 3)Place the hole in the base of the Stirrup over the footrest and wedge and push down until secure.
- 4)Take the strap on the Stirrup base, wrap around the leg to the footrest table leg, and attach to the velcro tab on the other side of the Stirrup.
- 5)Loosen the patient strap, have patient put their leg in the Stirrup and secure the patient strap over their leg. (optional)
- 6)Now the patient's can comfortably and securely rest their legs during their procedure or exam.
- 7)Reverse the order for **removal**, with care to avoid dropping the Stirrup onto your foot or floor- it is made of a durable product, but can be broken with misuse. It can remain attached to your table or be removed for cleaning, transport, or storage.
- 8)**Care of the Stirrup:** the entire unit is dish washer safe, and should be regularly cleaned and disinfected. Pads should be wiped down between uses. If Pads or Straps become soiled, they may be washed in a washing machine with mild detergent or hand washed. With regular use, Pads and Straps will become worn or torn. Replacements are available at davidscottco.com