



Cleaning Instructions: To ensure prolonged use of David Scott Company pads and straps please follow the below cleaning protocols.

1. Inspect pad before and after use to make sure there are no compromising defects
2. Use of disinfectant solution OTHER than those listed below may cause discoloration or deterioration of pad surface:
 - Coverage® Plus Germicidal Surface Wipes
Disinfecting/Deodorizing/Cleaning Wipes
 - Quaternary Ammonium Compound (Quats) With Ethanol Solvent
 - Quaternary Ammonium Compound (Quats) with Isopropyl Alcohol solvent
 - Quaternary Ammonium Compound (Quats) with IPA + 2-Butoxyethanol solvent
 - Quaternary Ammonium Compound (Quats) + Biguanide
 - Neutral Cleaners
3. The Use of H₂O₂ + PAA (Hydrogen Peroxide + Peracetic Acid) is strongly discouraged
4. **After using cleaner, wipe down with clean, damp, cloth with water.**
5. Always Follow manufacturer's instructions for concentrations and use of cleaning products.
6. Ensure no excess fluids remain on pads during and after cleaning
7. If recommended cleaning products are not available in your area, use a neutral disinfectant or disinfectant with sodium hypochlorite solution (.1% or 1000 parts per million available chlorine)



For special considerations such as large volume blood spills, contamination or exposure to high-risk tissue and fluids, pads may be cleaned with an intermediate to high level disinfection dilution rate of 10% or 1:10 dilution rate of household bleach as recommended by CFC and AORN guidelines. Please note this is not recommended for routine use and requires a damp cloth rinse post cleaning

Additional Precautions:

- Dilute all cleaners/disinfectants and germicides in accordance with manufacturer's instructions
- Allow adequate drying time before placing pad into service
- Some cleaning agents can leave a residue, it is recommended to rinse and wipe clean with a slightly damp cloth
- Avoid applying strong adhesive tape directly to mattress surface
- Use a soft sponge, cloth or wipe to clean; avoid hard bristle or abrasive brushes
- Never clean mattresses or pads with harsh cleaners or solvents that include moderate to high concentrations of alcohol
- Do not use cleaning agents designed for use on hard, non-porous or metal surfaces
- Do not allow any cleaning agent to remain in contact with fabric for a prolonged time or beyond cleaning agent manufacturer's recommendations
- Adequately rinse cleaning agents and disinfectants from fabric
- Final step should always be wiping down pad with a clean, damp cloth with water.

2 MASTER DRIVE FRANKLIN MA 02038
P: 508-875-3333 F: 508-875-3375